



USE INSTRUCTIONS

(α) Non Sterile Prosthetic Components

where α = FLY LINE, BUTTERFLY LINE, ZITAN LINE, SMILE LINE, PEARL LINE, DIAMOND LINE, OPAL LINE, CRYSTAL LINE, YALOS LINE, CARE LINE, WORLD LINE, SUN LINE, MOON LINE, LIFE LINE, WELL LINE, FAIR LINE, KING LINE, PRINCE LINE, MARK LINE

The non-sterile prosthetic components have been designed to be used for prosthetic equipment system ZivEco. For a full detail of the components, see the catalog of the manufacturer.

1- The non-sterile prosthetic components consist in:

- A) Castable abutment
- B) Clinched stump seed
- C) Transfer
- D) Similar or replication implant
- E) Instrumentation in general, except drills belonging to the rotary instruments.

2 -The packages are supplied sealed and integer; in case of damage to the packaging it's necessary to contact the seller or manufacturer for replacement.

OPERATING INSTRUCTIONS

Warning: This product has been designed and built exclusively for the provided conditions, any other use is considered improper and therefore dangerous.

The manufacturer cannot be held responsible for any damage caused by improper, erroneous and unreasonable; the product should be handled and used only by qualified staff, surgeons and dentists. In the case of defective products the manufacturer is exclusively responsible for the replacement of the same.

GENERAL

The twist of the prosthetical part must be conducted using the appropriate tool.

We recommend using the dynamometer ratchet exclusively for the final assembly of the prosthesis, not during the tests to avoid damages or breaks of the screws of the components. Use manual brucoles and always make sure the component and any accompanying screws are connected properly to prevent movement of lever, that increases the risk of fracture.

The reuse of products marked as single use can cause contamination of the product, infection to the patient and/or function of the device other than as provided.